



“Results matter.”

## VITROS® ECIQ

Immunodiagnostic  
System



Designed to help you  
ensure the result integrity  
& quality you need to  
confidently deliver the  
results that matter.

# A Versatile System that Provides Peace of Mind

## Advanced Security

- + The VITROS® ECiQ Immunodiagnostic System has proprietary Intellicheck® technology designed to significantly reduce analytical errors.
- + The IntelliReport provides real-time status on the quality of every patient result, for worry-free result reporting.

## Optimized Efficiency

- + Minimal maintenance, including no tubing cleaning or replacement, and no priming of reagents, even for start-up operations, means less operator intervention and true walk-away automation.
- + e-Connectivity® with Predictive Technology operations provide highly secure access to system data for immediate diagnostic and performance analysis, to maximize system uptime and increase staff productivity.

## Cost Effectiveness

- + Patented VITROS® Technologies are designed to deliver predictable financials and eliminate many of the hidden costs that can be associated with other immunoassay analyzers.
- + Standardized reagents use the same reagent packs across all VITROS® Systems

## Key Features

- + Enhanced Chemiluminescence technology for excellent sensitivity and broad dynamic ranges
- + MicroWell technology allows for small sample volumes with minimal waste
- + Intellicheck® technology for continuous process verification help ensure correct results the first time
- + Short sample, clot, and bubble detection offers greater security
- + Single use disposable tips eliminate carryover and cross contamination
- + No plumbing, drains, vents or deionized water eliminates many costs associated with alternative systems
- + Small footprint requires less space and overhead expenses
- + Up to 20 Integrated Reagent Packs loaded any time and stored on-board for continuous, random access operation
- + Low total costs with increased reportable result efficiency, fewer errors and interventions, and more effective staff utilization
- + Easy to train and cross-train operators
- + Up to 90 results per hour\*

\* Depending on test mix



Intellicheck®



MicroWell

AFP  
Anti-HAV IgM  
Anti-HAV Total  
Anti-HBc  
Anti-HBe  
Anti-HBs  
Anti-HBc IgM  
Anti-HCV  
Anti-HIV 1+2  
CA 125 11™  
CA15-3™  
CA19-9™  
CEA  
CK-MB  
CMV IgG°  
CMV IgM°  
Cortisol (serum & urine)  
Estradiol  
Ferritin  
Folate  
Free PSA°  
Free T3  
Free T4  
FSH  
HBeAg  
HBsAg  
HBsAg Confirmatory  
Intact PTH  
LH  
Myoglobin  
NT-proBNP  
NTx  
PSA  
Progesterone  
Prolactin  
RBC Folate

Rubella IgG  
Rubella IgM°  
Syphilis°  
Testosterone  
Total β-hCG II  
Total PSA II°  
Total T3  
Total T4  
Toxoplasma IgG°  
Toxoplasma IgM°  
TSH 3rd Gen  
T3 Uptake  
Troponin I ES  
Vitamin B12  
25-OH Vitamin D Total

HIV and Hepatitis co-developed with Novartis Vaccines and Diagnostics, Inc.

CA 125 II, CA 15-3, and CA 19-9 are trademarks of Fujirebi Diagnostics, Inc.

° Available outside the United States

Ortho Clinical Diagnostics

VITROS® System  
ECiQ  
Immunodiagnostic